CLINICAL IHT

24 chemin de l'Echut

31770 COLOMIERS

Tél : 05.62.74.08.44

Fax : 05.62.74.08.43

E-mail : clinical@clinical-iht.com

Le 15 juillet 2021

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EN

INDICATIONS:

Drenofast is a Redon drainage system with built-in vacuum designed to efficiently and continuously evacuate wound exudates post-operatively, with the purpose of facilitating the healing process and avoid complications, stimulating formation of granulation tissue and cell proliferation, reducing bacterial load and enhancing contraction of wound edges.

KONTRAINDICATIONS:

Catheters should not come into contact with tissue or organs that can be damaged by the vacuum system.

Drenofast should not be used in areas such as the brain, intestine, peritoneum, bone marrow or cavernous bone.

In case of use close to the dura mater and the possible risk of rupture: systematically check the condition of the dura mater and the compatibility of the desired depression. Do not use the Drenofast bottle without vacuum and do not puncture the bellows.

INSTRUCTIONS FOR USE:

- Open the container in sterile conditions.
- Remove the bottle from the protective bag.
- Remove the tube from the protective bag.
- Cut the sleeve at the place corresponding to the CH indicated for the catheter or adjust in case of screw connection.
- Connect the tube luer to the bottle, close the distal clamp and connect the catheter to the sleeve.
- Open the clamp on the bottle and the distal clamp. Suction will start.

CHANGE OF BOTTLE:

- Close the distal clamp.
- Close the bottle clamp.
- Disconnect the luer.
- Connect to the new bottle and open the bottle clamp.
- Open the distal clamp.

PRECAUTIONS:

Do not store at temperatures above 40°C.

Do not resterilise. The product must not be resterilised as this would affect the integrity of the materials, thereby favouring contamination and reducing their safety or efficacy. The accessories of the Drenofast product are designed to be used only with Redon type drainage systems with built-in vacuum.

It must be used by specialised health personnel. Prior to use, inspect the container and do not use the product if the container is open or damaged.

Prior to use, it is important to regularly check that the vacuum gauge is in the maximum position. If the bellows indicates minimum, change the bottle.

Avoid serous cavities and open wounds. Carefully remove and handle the needle cap and the needle to prevent possible cuts. Curved needles with double bevel tip present a greater risk of tearing the patient's tissue, consequently they should be used with care.

Record the implanted length of the catheter in the patient's clinical notes.

Secure the system connections.

The volume indicated on the bottle is approximate, maximum variation 5%.

During use, regularly check the correct functioning of the system. If the bellows indicates minimum, change the bottle.

If the vacuum bottle needs to be changed, first close the drain using the connection tube clamp to avoid a secretion reflux and then close the clamp of the bottle before disconnection.

Verify the total extraction of the catheter.

Destroy the device according to the healthcare centre's procedure for biological waste. This is a single-use device that is not designed to work as expected after the first use. Do not reuse, reprocess or resterilise. Reuse, reprocessing or resterilisation can compromise the structural integrity of the product, create risk of contamination and/or cause infection.

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